

**Accu-Chek Softclix**  
Nježan, praktičan  
i pouzdan.



## Accu-Chek Softclix

### Odabir dubine penetracije



Odaberite jednu od 11 postavki

### Umetanje lancete



Skinite kapicu



Umetnite lancetu



Odvrnite kapicu za sterilnost



Zamijenite kapicu

### Uzimanje uzorka krvi



Pritisnite gumb za pripremu



Pritisnite gumb za otpuštanje

### Zamjena lanceta



Skinite kapicu Izbacite iskorištenu lancetu



Umetnite novu lancetu



Odvrnite kapicu za sterilnost



Zamijenite kapicu



### Doživite što je sve moguće.

Da biste saznali više o proizvodima Accu-Chek i kako oni doprinose poboljšanju liječenja dijabetesa, posjetite [www.accu-chek.com](http://www.accu-chek.com).

ACCU-CHEK, CLIXMOTION, FASTCLIX, ACCU-CHEK MULTICLIX i SOFTCLIX žigovi su društva Roche. © 2015. Roche Diabetes Care

Svi drugi nazivi proizvoda i žigovi pripadaju odgovarajućim vlasnicima.

Roche Diabetes Care GmbH, Sandhofer Str. 116, D-68305 Mannheim